۲

CIRRUS[™] photo Scan Acquisition (Models 600 and 800) **Quick Reference Guide**

WARNING - GENERAL HAZARDS



(B Please read the CIRRUS photo user manual carefully before using this product.

C 1

CIRRUS CIRRUS photo Photo Select this mode for Select this mode Select this mode	 Select the desired imaging mode. Advanced tip: For patients who have difficulty fixating, take OCT and fundus images as
OCT imaging only for combined OCT for fundus imaging + fundus imaging only.	separate acquisitions (CIRRUS and Photo), and combine them using Review Image or Analyze
Align and fo cus the instrument using working distance dot	ts
CIRRUS photo File Tools Help	• Pull the instrument back from the patient and place patient's chin in the chinrest.
🔥 Patient 🙆 Acquire data 🚦 Review image I	Adjust table height and use eye level marks to adjust height of chinrest.
Tran, Lan 3. March 1933 CIRRUS CIRRUS photo Photo • OCT son Mediate cales 12/2021 •	• Use joystick to align instrument with patient pupil (use alignment circles for best positioning).
false correct	• Move instrument towards patient. The focusing aid and working distance dots (two white dots) will appear together with the fundus image.
	Bring the working distance dots roughly into focus.
	 If AutoFocus fails, use knob on side of instrument (bars should be aligned horizontally).
	 Bring the working distance dots into focus by turning and moving the joystick and align the working distance dots within the horizontal lines.
Acquire image(s)	
	 Instruct patient to blink once and to look directly at the fixation target Click joystick button to capture image(s). Advanced tip: To modify OCT scan pattern location, press arrow (•) on keypad before capturing images. Use controls on OCT fundus image to change the scan location and rotation/size (HD Raster only). Image: Use controls on OCT fundus image to change the scan location and rotation/size (HD Raster only).
	Grab and drag scan pattern to preferred scanGrab edge of scan pattern to resize the raster scan length and spacing. For 1 Line Raster, drag spacing to zero.Press on keypad or the Optimize button in display. Use the scroll wheel of your mouse to center the scan.
Review the captured scan, then save or try again.	
	• Evaluate the B-scan quality and the signal strength (should be 6 or greater).
Riceracy Luip 21. Dennier Pacific Status Sta	• Check that the OCT B-scans in the display windows appear centered and that no
	data is cut off or missing.
	Click Try Again to repeat scan acquisition, if necessary.
	Click Save to save the OCI scan (fundus images are automatically saved).



۲

۲





Carl Zeiss Meditec AG Goeschwitzer Str. 51-52 07745 Jena Germany
 Phone:
 +49 3641 220 333

 Fax:
 +49 3641 220 112

 Email:
 info.meditec@zeiss.com

 Internet:
 www.zeiss.com/med

000000-1795-718-KurzGA01-GB-280815 CIRRUS photo © Carl Zeiss Meditec AG, Jena Specifications subject to change

۲

۲



۲