



VISUPLAN® 500

The ZEISS VISUPLAN 500 is a non-contact tonometer to provide professional glaucoma care.

[Highlights](#)[Details](#)[Contact](#)

ZEISS VISUPLAN 500 enables you to provide easy glaucoma screening simply via a soft air puff – without contact and anaesthesia of the eye.

Easy handling

ZEISS VISUPLAN 500 is intuitive to handle with an automatic measuring process that is easily controlled via touchscreen.

A test air puff prepares your patient for the subsequent examination.

Fast and reliable examination

A stable headrest and LED fixation allow for a fast and reliable examination of your patient's intraocular pressure. The automatic positioning and measurement make sure your measurement results are repeatable and comparable. The results are directly displayed on the screen.

Benefits in glance:

- **Automatic positioning and measurement**
Thanks to the auto-tracking of the probe, it is very easy to find the right position and initiate the measurement.
- **Soft puff of air for more comfort**
The ZEISS VISUPLAN 500 makes the measurement of intraocular


pressure gentler and more comfortable for your patients. The non-contact measurement uses one or more puffs of air. Patients at risk of glaucoma can be quickly and reliably identified.

- **Stable headrest and reliable LED fixation**

The headrest easily moves, allowing for simple positioning of the patient. The integrated LEDs help patients intuitively fixate, enabling a fast and reliable measurement. Move the headrest to measure the second eye.

The key benefits for you

- You identify patients at risk quickly and reliably
- You receive accurate measuring results regardless of the user
- Your staff can take the measurement



Featuring measurement histogram to monitor the overall distribution.

ZEISS is one of the world's leading manufacturers of spectacle lenses, and is committed to delivering maximum precision and comfort. ZEISS designs and produces lenses, instruments and measurement systems, as well as retail concepts and technology services that continue to raise the bar in vision care.



Contact us to Get Started

ZEISS Customer Service

+44 (0)121 356 5595

Email



Explore

- Understanding Vision
- Health + Prevention
- Lifestyle + Fashion
- Driving + Mobility
- Sports + Leisure
- Work Life

Help me choose

- Distance + Reading Glasses
- Varifocal Lenses
- Sunglasses
- Working Glasses
- Sports Glasses
- Glasses for Children
- Lens Coatings
- Contact Lenses
- Lens Cleaning
- At the Optician

Services

- My Vision Profile
- Online Vision Check

For Eye Care Professionals

- Instruments + Technologies
- ZEISS Spectacle Lenses
- ZEISS Cleaning Solutions

