



## The Single Vision lens portfolio by ZEISS

Now with 4 performance levels.

### ZEISS offers a wide range of single vision products to serve all wearer needs.

From lower-price conventional single vision lenses to the highest precision with ZEISS Single Vision Individual lenses.


Visual quality, cosmetic appearance and wearing comfort, i. e. weight, must all be taken into account in the design of single vision lenses. Nowadays, higher refractive indices for both plastic and glass lenses and aspheric/atoroidal surface designs make it possible to produce flat, light lenses which provide even patients requiring high prescriptions with outstanding visual and wearing comfort.

Single vision lenses are lenses with only one dioptric power. They are primarily used for single power distance or reading spectacles. As single vision reading lenses often cannot provide presbyopes with acceptable visual quality at intermediate distances, they are now increasingly being replaced by so-called proximity lenses for more than the close range.

### Benefits of ZEISS single vision lenses

- ✓ Better sight at night: ZEISS single vision lenses are available with > **i.Scription® technology** for even greater contrasts and color intensity.
- ✓ Access to leading edge technology: e. g. hardest ever > **ZEISS coating DuraVision® Platinum** and self-tinting PhotoFusion lenses.
- ✓ Better vision and larger range of products; products – from conventional to totally individualised – serve all wearer needs.

## Products





**ZEISS Single Vision Individual**  
The customised single vision lens

[Learn more](#)



**ZEISS Single Vision Superb**  
Widen your view

[Learn more](#)



**ZEISS Single Vision AS**  
The Classic Plastic Lens with an Aspheric Surface Design

[Learn more](#)



**Single Vision Sph**  
The Classic Spectacle Lens

[Learn more](#)

## Depending on the consumer need, ZEISS offers different performance levels

### ZEISS Single Vision portfolio

**ZEISS  
Single Vision  
Sph**  
★★

**ZEISS  
Single Vision  
AS**  
★★★

**ZEISS  
Single Vision  
Superb**  
★★★★

**ZEISS  
Single Vision  
Individual**  
★★★★★

### Optimum visual performance

Incorporation of all relevant individual data measured by the optician or optometrist.

■

### Clearer vision over the entire lens

Advanced freeform technology allows more precision.

■

■

### Thinner & flatter lenses

Flatter curvature and reduced magnifying.

■

■

■

### High quality vision

Optimised surface allows punctual precise optical imagery.

■

■

■

■

### Better night vision

Available with i.Scription® technology.

■

■

■

■

### Single vision glass lenses are available from ZEISS in many different versions.

Depending on the power you need, lenses made of different materials are available to provide you with an optimum pair of glasses with the thinnest, lightest lenses possible.

## Related Products



### Coatings by ZEISS

High quality lenses require high quality treatments

[Learn more](#)



### ZEISS EnergizeMe Spectacle Lenses

The first spectacle lenses made for refreshing after contact lens time.

[Learn more](#)



### ZEISS DriveSafe Lenses

One pair of glasses for everyday use and safer driving.

[Learn more](#)

ZEISS is one of the world's leading manufacturers of spectacle lenses, and is committed to delivering maximum precision and comfort. ZEISS designs and produces lenses, instruments and measurement systems, as well as retail concepts and technology services that continue to raise the bar in vision care.



### Contact us to Get Started

ZEISS Customer Service

[+44 \(0\)121 356 5595](#)

[Email](#)





---

**Explore**

- Understanding Vision
- Health + Prevention
- Lifestyle + Fashion
- Driving + Mobility
- Sports + Leisure
- Work Life

**Help me choose**

- Distance + Reading Glasses
- Varifocal Lenses
- Sunglasses
- Working Glasses
- Sports Glasses
- Glasses for Children
- Lens Coatings
- Contact Lenses
- Lens Cleaning
- At the Optician

**Services**

- My Vision Profile
- Online Vision Check

**For Eye Care Professionals**

- Instruments + Technologies
- ZEISS Spectacle Lenses
- ZEISS Cleaning Solutions